

Product info

Duren Medi-Line-PU-med

producer : Eurofoam GmbH, Greinerstraße 70, A-4550 Kremsmünster

Tel: +43 7583 90500 259, Fax: +43 7583 90500 389, e-Mail: foamforcare@eurofoam.at

1. composite

Duren Medi-Line-PU-med parts are made up of a top quality Polyurethan-foam, with a liquid refusing and almost airproof PU-coating.

2. purpose

Duren Medi-Line-PU-med are exclusive parts for standby, immobilisation, for bearing or fixation from patients (e.g.. half roll, three-quarter roll, roll, cube, bearing wedge, head support, rings, discs, burden, leg support). Use Duren Medi-Line-cell only for patients with entire skin.

3. classification

Duren Medi-Line parts are medicine products which rank among **class 1 (not steril, not invasiv)**.

4. features

- breathable, refuse liquid, single- or multilayer PU-skin
- supporting, high elastic bearing capacity
- easy to clean and hygienic, sturdy, formfitting adjustment
- not antistatic and not electroconductive
- free of latex

5. application- and cleaning advice

- using temperature +10 to +50 °C
- autoclavable to 135 °C saturated steam, surface disinfection with commercial disinfectant (notice the application at the particular product specification)
- enough cleaning after application
- maximal 15 passes at the steam autoclave
- storage of the parts only in dry, light proof and clean ambit

6. possible dangers if using Duren Medi-Line-PU-med :

- If diverting from its intended use or by incorrect, faulty application it comes to a loss of characteristics.
- Persons who are allergic or sensible to polyurethane should avoid the contact with Duren Medi-Line-PU-med parts.
- The influence of ignition sources or heat sources creates fire danger.
- Damage of the coating through the application, cleaning and disinfection; through sharp-edged objects
- Danger of drown by the almost airproof surface of the parts
- The almost liquid tight surface can't absorb perspiration.
- Decubitus danger by supporting the patient too long.
- Discoloration, yellowing and brittleness by influence of sunlight or UV-light.
- Lots of germs by unclean storage, insufficient cleaning or disinfection
- Damage by common use or application of solvent for cleaning.
- Reddened skin is possible
- Perspiration is possible

7. durability and disposal

Duren Medi-Line parts are by correct storage durable for 1 year.

Please dispose not contaminated articles with commercial waste. For contaminated articles it's the best to destroy them in waste combustion